BraceTek[™] Knee ROM and Back LSO

Provides Stability alongside Cold, Heat and Compression Therapy

ThermoTek's BraceTek options are made with quality materials to enhance comfort and support. Options are available for Knee ROM and Back LSO. The braces are easy to adjust and apply and are available with removable thermal/compression wraps for use with our VascuTherm devices.

Knee ROM Features

- Adjustable posterior hook and looplock closure (5" adjustable), available in 13"
- Comfortable large popliteal opening
- Polycentric flip-top hinge with easy set flexion/extension at 0°, 15°, 30°, 60° and 90°
- Metal double uprights for medial/lateral knee control
- Hinges can be inverted for future wear as a +rehab brace
- Upper and lower circumferential straps with elastic segments
- Indicated For:
 - + Preventive Care, Arthritis
 - + Post-Op management
 - + MCL/LCL instability
 - + PCL/ACL Grade I & II instability

Size	Knee Circumference	6" above Knee
XS	12-13"	16-19"
S	13-14"	17-20"
М	14-15"	18-21"
L	15-16"	19-23"
XL	16-18"	20-24"
2XL	18-20"	22-26"
3XL	22-24"	24-28"

ThermoTek

Backed by 30 years of experience with thermoelectric and precision engineering, ThermoTek offers advanced solutions that support post-surgical recovery and scientific thermal management. We pride ourselves on delivering high-quality dependable products that directly address our customers' needs.

Back LSO Features

- Provides coverage from sacral spine through T-9 vertebrae
- Sagittal control with rigid anterior and posterior panels.
- Lumbar pad and overlapping double pull strap design provides lumbar and abdominal support.
- Provides consistent pain relief.
- One-Size-Fits-All, adjustable with extender panels.
- Indicated For:
 - + Post-surgical Lumbar Stabilization
 - Herniated and Bulging Discs
 - + Lower and Mid-Back Pain
 - + Spinal Stenosis
 - + Spondylosis and Spondylolisthesis
 - + Acute/Chronic Strains
- **Review contraindications before use**

NOTE: Federal law restricts this device to use by or on the order of a licensed medical practitioner. Follow the prescribed instructions of your medical practitioner for therapy mode, treatment area, duration and frequency of treatment.









Think Customer Enable Change

Drive simplified solutions.

Challenge the status quo.

One team, one direction.

Alignment

Own it, do it right.

Empowerment Sustainability Accelerate design

and development

phone: 1-877-242-3232

www.thermotekusa.com

OPIFMTGBKRBL_A 04/24